


at last . . . a 400 tube fraction collector!

the new fractomette 400

*Collect the largest number of samples
in the least amount of space.*

If you wish to collect a large number of fractions from a fast column or wish to collect 200 fractions each from two columns simultaneously, the Buchler Fractomette 400 is ideal for you. Collection by Drop, Time or Volume in 18 x 150 MM test tubes. Equipped with automatic Drop-Stop preventing stray drops or loss of samples due to power failure. Sturdy aluminum frame for column mounting. This latest addition to Buchler's family of fraction collectors costs only \$1375.



MANUFACTURERS' LITERATURE

Spectroscopy Glassware. Bulletin No. 5006-1, 28 pages, presents sample tubes, coaxial cells, and other equipment and special devices for NMR and EPR spectroscopy. Wilmad Glass Co., Inc., Route 40 and Oak Rd., Buena, N. J. 08310 **601**

Variable Pathlength Cell. When and how to use variable pathlength cells in solvent compensation and differential analysis is discussed in a new data sheet. Barnes Engineering Co., 30 Commerce Rd., Stamford, Conn. 06902 **602**

Helium Ionization Detector. Bulletin 4084, two pages, describes a helium ionization detector for trace analysis of all inorganic as well as organic gases with a sensitivity of 1 to 2 ppm. Beckman Instruments, Inc., 2500 Harbor Blvd., Fullerton, Calif. 92634 **603**

Chemicals. Aliphatic monocarboxylic, dicarboxylic, aromatic, and fatty acids are included in this list of organic chemicals. Chemical & Radioisotopes Div., International Chemical & Nuclear Corp., 801 N. Lake St., Burbank, Calif. 91502 **604**

Mass Spectrometers and High Vacuum Equipment. Residual gas analyzers, mass spectrometers, traps and accessories, and UHV systems are presented in a 16-page catalog. Aero Vac Corp., P. O. Box 448, Troy, N. Y. 12181 **605**

Dissolved Oxygen Analyzer. Model 735 dissolved oxygen analyzer and its advanced polarographic sensor are described in Bulletin 4092, four pages. Beckman Instruments, Inc., 2500 Harbor Blvd., Fullerton, Calif. 92634 **606**

Ge(Li) Detectors. "The Selection and Use of Ge(Li) Detectors" is a 16-page manual designed as a guide for the researcher. Nuclear Diodes, Inc., P. O. Box 135, Prairie View, Ill. 60069 **607**

Fractionation of Collagen Subunits. Data Sheet No. 3 describes the isolation of collagen $\alpha 1$, $\alpha 2$, $\beta 11$, and $\beta 12$ from heat denatured, salt-extracted and acid-extracted rat skin collagens by using Whatman CM32. Reeve Angel, 9 Bridewell Place, Clifton, N. J. 07014 **608**

Water Quality Monitor. The features of Model 9500, which can monitor up to nine parameters, are presented in Bulletin 4083, eight pages. Beckman Instruments, Inc., 2500 Harbor Blvd., Fullerton, Calif. 92634 **609**