

- July 17 to 21—Industrial Use of the Polarizing Microscope.** London, England. McCrone Research. £55. Contact (7)
- July 17 to 21—Forensic Document Microscopy.** Georgetown University. J. M. England. Contact (5)
- July 20 to 21—Color Technology for Management.** Rensselaer Polytechnic Institute. Fred W. Billmeyer, Jr. Contact (9)
- July 24 to 28—Applied Molecular Spectroscopy: Infrared, Raman, Ultraviolet.** Arizona State University. \$200. Contact: Jacob Fuchs, Modern Industrial Spectroscopy, Arizona State University, Tempe, AZ
- July 24 to 28—Advances in Color Technology.** Rensselaer Polytechnic Institute. Fred W. Billmeyer, Jr. Contact (9)
- July 24 to 28—Forensic Analysis of Writing Inks.** Georgetown University. E. P. Cofield, Jr. (Scripto, Inc.); R. W. McKinney (W. R. Grace & Co.). Contact (5)
- July 24 to 28—Forensic Microscopy of Human and Animal Hairs and Fibers.** Georgetown University. Marie Jones (Bristol-Myers Co.). Contact (5)
- July 24 to Aug. 4—Forensic Optical Microscopy.** Georgetown University. Mary L. Willard (Pennsylvania State University). Contact (5)
- July 27 to 28—Thin-Layer Chromatography Workshop.** Madison, WI. \$125. Contact: James E. Nicholls, Dept. of Engineering, University of Wisconsin Extension, 432 Lake St., Madison, WI. 608-262-2061
- July 31 to Aug. 4—Two Separate Courses: Laboratory Management; Gas Chromatography.** Philadelphia, PA. Sadtler Research. \$275 each. Contact (3)
- July 31 to Aug. 4—Forensic Analysis of Ethyl Alcohol.** Georgetown University. Dr. H. C. Freimuth (Office of the Chief Medical Examiner, Baltimore, MD). Contact (5)
- July 31 to Aug. 4—Forensic Microscopy of Synthetic and Cellulosic Fibers.** Georgetown University. R. Scott (Du Pont). Contact (5)
- July 31 to Aug. 4—Infrared Spectroscopy: Applications.** Bowdoin College, ME/MIT. \$300 with room and board. Contact: Dana W. Mayo, Bowdoin College, Brunswick, ME 04011. 207-725-8731, X601
- Aug. 7 to 18—Modern Industrial Spectroscopy.** Tempe, AZ. \$300. Contact: Jacob Fuchs, Director, Modern Industrial Spectroscopy, Chemistry Dept., Arizona State University, Tempe, AZ 85281
- Aug. 7 to 18—Forensic Microscopic Analysis of Inorganic Ions.** Georgetown University. Elizabeth Anne Whitman (Naval Ordnance Station, Indian Head, MD). Contact (5)
- Aug. 14 to 18—Fisk Institute: Two Separate Courses: Basic Infrared Spectroscopy; Gas-Liquid Chromatography.** Nashville, TN. \$200 each. Contact (10)
- Aug. 21 to 25—Fisk Institute: Two Separate Courses: Practical Laser Raman Spectroscopy; Interpretation of Infrared and Raman Spectra.** Nashville, TN. \$200 each. Contact (10)
- Sept. 4 to 8—Reflected Light Microscopy.** London, England. McCrone Research. £55. Contact (7)
- Sept. 11 to 13—Analytical Chemistry Workshop.** Ottawa, ONT. The Chemical Institute of Canada. A. H. Gillieson, S. Nargowalla, R. F. Pottie, A. G. Plant, J. A. Page. \$25. Contact: M. E. Plummer, Department of Chemical, Biochemical, and Medical Technologies, Algonquin College, 200 Lees Ave., Ottawa, ONT, Canada
- Sept. 18 to 22—Color and the Behavior of Colorants.** Attleboro, MA. Kollmorgen. \$250. Contact (6)
- Oct. 4 to 6—Interpretation of Mass Spectra.** Holiday Inn South, Indianapolis, IN. F. W. McLafferty. \$135. Contact: R. M. Brooker, Center for Industrial Relations and Continuing Education, Indiana Central College, 4001 Otterbein Ave., Indianapolis, IN 46227

- (1) Education Dept., American Chemical Society, 1155 Sixteenth St., N.W., Washington, DC 20036. Late registrants, please call 202-737-3337, X258
- (2) Camag, Inc., 2855 S. 163rd St., New Berlin, WI 53151. 414-786-9550
- (3) Sadtler Educational Div., Sadtler Research Laboratories, Inc., 3316 Spring Garden St., Philadelphia, PA 19104. 214-382-7800
- (4) Washington University, Box 1048, St. Louis, MO 63130. 314-863-0100
- (5) Joseph M. English, Forensic Sciences Center, 9 Reiss Science Bldg., Georgetown University, Washington, DC 20007
- (6) Ruth M. Johnston, Kollmorgen, Color Systems Div., 67 Mechanic St., Attleboro, MA 02703. 617-222-3880
- (7) Mrs. Miriam L. Fallert, McCrone Research Institute, 451 E. 31st St., Chicago, IL 60616. 312-842-7105
- (8) Martin Hughes, Varian Assoc. of Canada, Ltd., 6358 Viscount Rd., Malton, ONT, Canada
- (9) Office of Continuing Studies, Color Technology Program, Rensselaer Polytechnic Institute, Troy, NY 12181. 518-270-6531
- (10) Enrique Silberman, Acting Director, Fisk Institute, Box 8, Fisk University, Nashville, TN 37203

**entific Instrument Co.,** 69-57 Juniper Blvd. South, Middle Village, NY 11379. This company will operate as a division under Spectrum's Laboratory Products Division. Post makes countercurrent distribution instruments, counter-dou-

ble-current distribution instruments, and countercurrent thin-film dialysis equipment.

**Rigaku Denki Co., Ltd.,** has announced the formation of **Rigaku/USA.** The

## From Manostat's line of quality products...

REPEATING PIPETTERS FOR ACIDS AND SOLVENTS

### Minipet®

Adjustable syringe type dispenser featuring two-way check valve. Just push to dispense and release to fill from remote containers. Teflon® plunger and glass barrel assures non-sticking smooth operation. Autoclavable. Accuracy.  $\pm 0.5\%$  repetitive.

Sizes. 1, 2, 10 and 30 ml.

### Varipet®

Variable volume transfer pipet with similar precise, chemically inert and non-sticking construction features as Minipet®.

Availability. 0.5 & 1 ml (0.01 ml divisions), 2 & 4 ml (0.05 ml divisions), 10 ml (0.2 ml divisions) and 30 ml (1 ml divisions). Metal luer lock, glass luer or Teflon tips.

ASK FOR DETAILED BULLETIN

Manostat products are available from your laboratory supply dealers


**Manostat**  
20 N. Moore St., New York, N.Y. 10013

FOR LITERATURE CIRCLE NO. 135  
FOR DEMONSTRATION CIRCLE NO. 136

## From Manostat's line of quality products...

ACCURATE AUTOMATIC DILUTER FOR RATIOS UP TO 1000:1

### Dilutrol™



Fingertip controlled diluter provides precise pickup of predetermined amount of sample. Unique airlock between sample and diluent prevents cross contamination. Rinsing action of diluent cleans tip for next cycle. Liquids touch only Teflon® and glass. Dispensing head removable for operation up to 4 feet away.

Accuracy. Sample repeatability within  $\pm 0.3\%$ . Diluent delivery to  $\pm 0.5\%$ .

Sizes. 0.5, 1 and 4 ml adjustable sample syringes each with 10 ml diluent syringe.

ASK FOR DETAILED BULLETIN

Manostat products are available from your laboratory supply dealers

**Manostat**  
20 N. Moore St., New York, N.Y. 10013

FOR LITERATURE CIRCLE NO. 137  
FOR DEMONSTRATION CIRCLE NO. 138