

FUTURE ARTICLES

Listed below are papers to appear in forthcoming issues of *Analytical Chemistry*

A New Concept in Precision Photometric Analysis Using a Radioisotopic Light Source

H. H. Ross

Analysis of Alkylene Oxide Polymers by Nuclear Magnetic Resonance Spectrometry and by Gas Liquid Chromatography

A. Mathias and N. Mellor

Automatic Continuous Determination of Small Quantities of Gaseous Ammonia

N. S. Segal and R. Wodley-Smith

Determination of Coal in Formalin-Fixed Pneumoconiotic Lungs: A New Method of Tissue Digestion Using Glacial Acetic Acid

I. Bergman

Determination of Silver in Minerals and Ores by Neutron Activation Analysis and High Resolution Gamma Spectrometry

G. L. Schroeder, R. D. Evans, and R. C. Ragaini

Critical Study of Three Methods for the Determination of Oxygen in Thorium and Yttrium Metals

V. A. Fassel, W. E. Dallmann, and C. C. Hill

Comparative Study of the Kinetics of the Benzidine Rearrangement by Four Electrochemical Techniques

D. M. Oglesby, J. D. Johnson, and C. N. Reilly

The Separation of Pb(II) From Bi(III), Tl (III), Cd(II), Hg(II), Au(III), Pt(IV), Pd(II), and Other Elements by Anion Exchange Chromatography

F. W. E. Strelow and F. V. Toerien

Titrimetric and Equilibrium Studies Using Indicators Related to Nile Blue A

M. M. Davis and H. B. Hetzer

Relationship between Resonance Line Profile and Absorbance in Atomic Absorption Spectrometry

K. Yasuda

Syringe Ram for a Rapid-Freeze Sampling Instrument

R. E. Hansen and H. Beinert

Determination of Emissivities with a Differential Scanning Calorimeter

R. N. Rogers and E. D. Morris, Jr.


Estimating Activation Energies with a Differential Scanning Calorimeter

R. N. Rogers and E. D. Morris, Jr.

For more accurate Absorption Spectrophotometry SHIMADZU PHOTOELECTRIC SPECTROPHOTOMETER

Shimadzu

Shimadzu QV-50 Photoelectric Spectrophotometer features a Littrow type Monochromator utilizing a fused silica 33 degree prism. Evacuated by a vacuum pump, it is capable of measurement in the extreme ultraviolet, down to 183 m μ . Constant temperature device is built in the cell compartment. More efficient for biological research. Especially flow dichroism measuring set and circular dichroism measuring set, both developed by Shimadzu, and other various accessory sets—Absorption Tube Accessory Set, Flame Photometry Accessory Set, Reflection Accessory Set, and Fluorometry, Nephelometry Accessory Set—are useful optional equipment to the QV-50.



SHIMADZU SEISAKUSHO LTD.

FOREIGN TRADE DEPT.

TOKYO OFFICE: Kanda-Mitoshicho, Chiyoda-ku, Tokyo

(Cable Add.: SHIMADZU TOKYO)

KYOTO OFFICE: Nishinokyo-Kuwabaracho, Nakagyo-ku, Kyoto

(Cable Add.: SHIMADZU KYOTO)

Sole Agents in U.S.A.

ATAKA NEW YORK INC.

633, Third Avenue, New York 17, N.Y.

Tel. No. OXford 7-7480

Circle No. 225 on Readers' Service Card